



Legionella Chain of Custody

LA Testing Order Number / Lab Use Only

PHONE:
EMAIL:

If Bill-To is the same as Report-To leave this section blank. Third-party billing requires written authorization.

Customer Information	Customer ID:		Billing Information	Billing ID:	
	Company Name:			Company Name:	
	Contact Name:			Billing Contact:	
	Street Address:			Street Address:	
	City, State, Zip:	Country:		City, State, Zip:	Country:
	Phone:			Phone:	
	Email(s) for Report:			Email(s) for Invoice:	

Project Information			
Project Name/No:		Purchase Order:	
LA Testing LIMS Project ID: <small>(If applicable, LA Testing will provide)</small>	US State where samples collected:	State of Connecticut (CT) must select project location: Commercial (Taxable) Residential (Non-Taxable)	
Sampled By Name:	Sampled By Signature:	No. of Samples in Shipment	

Sterile, Sodium Thiosulfate Preserved Bottled Used	Biocide Used in Source (If Checked, Please Specify Type)
Bottles Chilled on Frozen Ice Packs	Public Water Supply Samples Note: All results may automatically be reported to DOH if required by State.

Legionella Analysis (Performed by ISO 11731:2017 Culture Method) - TAT's are subject to methodology requirements.

M341: *Legionella pneumophila* confirmation and enumeration. Individual serotyping not included. Other species reported as *Legionella sp.*

M342: *Legionella pneumophila* confirmation and enumeration of individual serotypes 1-14. Other species reported as *Legionella sp.*

M343: *Legionella* Species confirmation and individual enumeration of 11 Legionella species (L. anisa, L. bozmaniai, L. dumoffii, L. feeleeii, L. gormanii, L. jordanis, L. longbeachae, L. maceachernii, L. micdadei, L. pneumophila, and L. sainthelensi) and their individual serotypes. Other species can be requested at an additional charge.

M344: ISO filter direct for water expected to have low bacterial concentrations (bottled water water for dialysis, distilled water.) 11 *Legionella* species (above) and their respective serotypes, identified and enumerated individually. Other species reported as *Legionella sp.*

M345: Identification and enumeration of *L. pneumophila* serotype 1, *L. pneumophila* serotypes 2-14, and *Legionella sp.*

M214: Pure Culture Preparation and Storage

M215: Pure Culture Preparation and Shipping

Heterotrophic Plate Count (HPC) Analysis (Check one below if requested) Test Code M015 Heterotrophic Plate Count (HPC) Standard Method 9215

Plate count agar, 35°C, 48 hour incubation R2A agar, 25°C, 5-7 day incubation

Sample #	Sample Location	Sample Type	Potable / Non-Potable	Test Code	Volume	Date / Time Collected	Temperature (Lab Use Only)
Example: Sample 1	Kitchen	Water	Potable	M343	1,000 ml	1/1/2021 3:30pm	

Special Instructions and/or Regulatory Requirements (Sample Specifications, Processing Methods, Limits of Detection, etc.)

Method of Shipment:		Sample Condition Upon Receipt:		Received on Ice? Check if Yes:	
Relinquished by:	Date/Time:	Received by:	Date/Time:	Received by:	Date/Time:
Relinquished by:	Date/Time:	Received by:	Date/Time:	Received by:	Date/Time:

Controlled Document - COC-26-LAT Legionella LA Testing R6 11/07/2023

AGREE TO ELECTRONIC SIGNATURE (By checking, I consent to signing this Chain of Custody document by electronic signature.)

EMSL Analytical, Inc.'s (DBA LA Testing) Laboratory Terms and Conditions are incorporated into this Chain of Custody by reference in their entirety. Submission of samples to LA Testing constitutes acceptance and acknowledgment of all terms and conditions by Customer.

